



Product Specifications

Defibrillator AED Awareness Bundle Pack

- Step by step instructions on defibrillation and the use of AED unit
- Increase visibility of the AED unit and raise awareness of the trained personnel

The Defibrillator AED Awareness Bundle Pack is ideal for increasing the awareness of the Defibrillator Heart Starter unit, its location and the trained personnel. The Defibrillator AED Awareness Bundle Pack also provides step by step instructions on defibrillation, CPR and how to use an AED unit for both trained an untrained personnel. The pack contains the following items:

- AED for Untrained Personnel Poster
- AED Defibrillation and CPR Poster
- 2 x Automated External Defibrillator Trained Personnel Sign
- 2 x Automated External Defibrillator Sign
- 2 x Automated External Defibrillator Nearest AED Sign

AED Defibrillation & CPR

1. AED Defibrillation

What is an AED?
An AED is a portable, computerized medical device that can be used to deliver an electric shock to a person's chest to restart a stopped heart.

How does it work?
The AED will analyze the heart's rhythm and determine if a shock is needed. If a shock is needed, the AED will deliver a shock to the heart to restart it.

When should you use an AED?
An AED should be used on anyone who is unconscious and not breathing normally.

2. Danger

Always use caution:
• Do not touch the person while the AED is analyzing the heart or delivering a shock.
• Do not touch the person while the AED is delivering a shock.

Always use caution:
• Do not touch the person while the AED is analyzing the heart or delivering a shock.
• Do not touch the person while the AED is delivering a shock.

3. Response

Check response:
• Tap and shout. If there is no response, shout for help.
• Check for breathing. If there is no breathing, shout for help.
• Check for pulse. If there is no pulse, shout for help.

Call for help:
• Call 999 or 112.
• Tell the operator you need an ambulance.
• Give the operator your location.
• Give the operator the name of the person who is unconscious.

4. Airway & Breathing

Check for airway:
• Tilt the head back.
• Lift the chin.

Check for breathing:
• Look for chest rise.
• Listen for breath sounds.
• Feel for airway resistance.

5. Getting Help

Call 999 or 112:
• Call 999 or 112.
• Tell the operator you need an ambulance.
• Give the operator your location.
• Give the operator the name of the person who is unconscious.

Call the Emergency Medical Services:
• Call the Emergency Medical Services.
• Tell the operator you need an ambulance.
• Give the operator your location.
• Give the operator the name of the person who is unconscious.

6. CPR

Check for response:
• Tap and shout. If there is no response, shout for help.
• Check for breathing. If there is no breathing, shout for help.
• Check for pulse. If there is no pulse, shout for help.

Call for help:
• Call 999 or 112.
• Tell the operator you need an ambulance.
• Give the operator your location.
• Give the operator the name of the person who is unconscious.

7. Using the AED

Check for response:
• Tap and shout. If there is no response, shout for help.
• Check for breathing. If there is no breathing, shout for help.
• Check for pulse. If there is no pulse, shout for help.

Call for help:
• Call 999 or 112.
• Tell the operator you need an ambulance.
• Give the operator your location.
• Give the operator the name of the person who is unconscious.

AED Automated External Defibrillators For Untrained Personnel

1. What is an AED?

What is an AED?
An AED is a portable, computerized medical device that can be used to deliver an electric shock to a person's chest to restart a stopped heart.

How does it work?
The AED will analyze the heart's rhythm and determine if a shock is needed. If a shock is needed, the AED will deliver a shock to the heart to restart it.

When should you use an AED?
An AED should be used on anyone who is unconscious and not breathing normally.

2. Who can use an AED?

Who can use an AED?
Anyone can use an AED. It is designed to be used by untrained personnel.

How to use an AED:
• Turn on the AED.
• Attach the pads to the person's chest.
• Follow the voice prompts.

3. How do AEDs work?

How do AEDs work?
The AED will analyze the heart's rhythm and determine if a shock is needed. If a shock is needed, the AED will deliver a shock to the heart to restart it.

When should you use an AED?
An AED should be used on anyone who is unconscious and not breathing normally.

4. When should an AED be used?

When should an AED be used?
An AED should be used on anyone who is unconscious and not breathing normally.

When should you use an AED?
An AED should be used on anyone who is unconscious and not breathing normally.

5. How to use an AED

How to use an AED:
• Turn on the AED.
• Attach the pads to the person's chest.
• Follow the voice prompts.

When should you use an AED?
An AED should be used on anyone who is unconscious and not breathing normally.



Defibrillator Heart Restarter



Defibrillator Heart Restarter

Trained Personnel Ext.



Defibrillator Heart Restarter

Nearest Defibrillator Heart Restarter is situated at

Person in charge