

FIRST



Product Specifications

Evolution LifeVac Complete First Aid Point, British Standard Compliant

- **Life-Saving Device:** Combines standard burn, biohazard and first aid kits with an advanced anti-choking device.
- **Regulated Medical Device** – CE marked, UKCA marked, MHRA regulated and MDR compliant
- **Includes a British Standard First Aid Kit** – A comprehensive and legally compliant First Aid Point for managing a wide range of workplace accidents and emergencies.
- **Addresses A Known Risk** - Provides a life-saving solution for choking emergencies that traditional kits don't cover.
- **Perfect for workplaces, schools, care homes, and public settings** seeking full emergency readiness.

The **Evolution LifeVac Complete First Aid Point** is a comprehensive first aid solution designed for busy workplace environments. It provides a centralised and easily accessible station to support effective response across a wide range of potential incidents. **Key Benefits**

- **Includes first aid kits to cover multiple scenarios** – this first aid point has the supplies to provide general first aid and burns first aid, deal with biohazards and remove blockages in a casualty's airway.
- **Includes a proven life-saving device** – the LifeVac anti-choking device has been proven to save lives in over 5,000 cases.
- **Interlocking, quick release brackets** – the brackets can be arranged in multiple configurations allowing the point to be constructed in a way that best fits the space available.
- **Compliant with British Standard 8599-1** – the general first aid kit that is included in the Evolution LifeVac Complete First Aid Point contains all the supplies that are required under British Standard 8599-1.
- **Photoluminescent headers and bracket labels** – improve visibility and readiness by ensuring first aid points remain easy to locate in low-light conditions and by helping staff quickly identify and replace any missing kits.
- **Regulated Medical Device** – CE marked, UKCA marked, MHRA regulated and MDR compliant
- **Supplied in Evolution cases** – the durable Evolution cases ensure that supplies are kept protected from damp and dust.
- **Enables swift intervention during critical incidents** – the LifeVac allows rapid action in choking emergencies, while the first aid kit provides essential supplies for cuts, burns, sprains, and other common injuries, helping to minimise complications.
- **Promotes employee confidence and peace of mind** – knowing that both everyday first aid and life-saving choking devices are readily available helps staff feel secure, empowered, and prepared to act in an emergency.

This First Aid Point Contains: BS 8599-1 First Aid Kit

- 1 x First aid guidance leaflet
- 4 x HypaCover First aid dressings, 12 x 12cm
- 3 x HypaCover First aid dressings, 18 x 18cm
- 3 x HypaBand Triangular bandages
- 3 x HypaCover Eye dressings
- 60 x HypaPlast Washproof plasters
- 30 x HypaClean Sterile wipes
- 2 x HypaPlast Microporous tape, 2.5cm x 5m
- 9 x HypaTouch Nitrile gloves (pair)
- 3 x HypaCover Finger dressings
- 1 x HypaGuard Face shield
- 2 x HypaGuard Foil blanket
- 2 x HypaSoothe Burn dressings, 10x10cm
- 1 x Clothing cutters
- 2 x HypaBand Conforming bandages, 7.5cmx4m

LifeVac Kit

- 1 x LifeVac User Guide Booklet
- 1 x LifeVac 'Instructions for Use' Leaflet
- 1 x LifeVac Unit.
- 1 x Paediatric Mask.
- 1 x Small Child/Adult Mask.
- 2 x Medium Adult Mask.
- 1 x Practice Mask.
- 1 x Large Adult Mask.
- 1 x A3 BLS Poster
- 1 x Window/Door Sticker

- 1 x Feedback Form

Burns Kit

- 1 x First aid guidance leaflet
- 2 x HypaSoothe Burn dressings, 10x10cm
- 3 x HypaSoothe Burn gel sachets, 3.5g
- 12 x HypaBand Safety pins
- 2 x HypaBand Conforming bandages, 7.5cmx4m
- 2 x HypaTouch Disposable Gloves (Pair)

Body Fluid Disposal Kit

- 2 x 10g HypaClean Absorbent powder
- 2 x 50ml HypaClean Disinfectant cleaner spray
- 2 x Disinfectant wipes
- 2 x HypaTouch Disposable gloves (pair)
- 4 x Non-woven cloths
- 2 x Scraper/scoop
- 2 x Polythene aprons
- 2 x Biohazard disposal bags

Ideal for:

- Care homes and assisted living facilities
- Healthcare clinics and rehabilitation centres
- Schools and childcare centres
- Workplaces with at-risk employees
- Community centres and public venues
- Factories / manufacturing plants
- Warehouses or distribution centres
- Construction sites
- Large offices or multi-storey buildings





First Aid Contents

First Aid Kit Contents

Items	Qty
First Aid Guidance Leaflet	1
HypaCover First Aid Dressings, 12cm x 12cm	4
HypaCover First Aid Dressings, 18cm x 18cm	3
HypaBand Triangular Bandages	3
HypaCover Eye Dressings	3
HypaPlast Washproof Plasters	60
HypaCover Sterile Wipes	30
HypaCover Microporous Tapes, 2.5cm	2
HypaTouch Nitrile Gloves (Pair)	9
HypaCover Finger Dressings	3
HypaGuard Face Shield	1
HypaGuard Foil Blankets	2
HypaSoothe Burn Dressings, 10cm x 10cm	2
Clothing Cutters	1
HypaBand Conforming Bandages, 7.5cm	2

Burns Kit Contents

Items	Qty
First Aid Guidance Leaflet	1
HypaSoothe Burn Dressings, 10x10cm	2
HypaSoothe Burn Gel Sachets, 6g	3
HypaBand Safety Pins	12
HypaBand Conforming Bandages, 7.5cm	2
HypaTouch Disposable Gloves (Pair)	2

LifeVac® Choking Rescue Contents

Items	Qty
Paediatric Mask	1
Small Child/Adult Mask	1
Medium Adult Mask	1
Large Adult Mask	1
LifeVac® Anti-Choking Device	1
Practice Mask (Medium Adult)	1
BLS Poster, A3	1
Window/Door Sticker	1
Instructions For Use Leaflet	1
LifeVac® Instruction Booklet	1

Body Fluid Disposal Contents

Items	Qty
HypaClean Absorbent Powders, 10g	2
HypaClean Disinfectant Cleaner Spray, 50ml	2
Disinfectant Wipes	2
HypaTouch Disposable Gloves (Pair)	2
Non-Woven Cloths	4
Scrapers & Scoops	2
Polythene Aprons	2
Biohazard Disposal Bags	2

INSTRUCTIONS FOR USE IN THE CASE OF CHOKING
If a casualty is a victim of choking

- 1 Encourage the victim to cough.
- 2 If the cough is not effective, call Emergency Services and give 5 Back Slaps.
- 3 If obstruction has not been cleared, follow by 5 Abdominal Thrusts.



Instructions For Use

1. Encourage the victim to cough. If not effective call the emergency services and follow by 5 back slaps and 5 abdominal thrusts. Do not apply LIFEVAC until obstruction has been cleared.

2. If the cough is not effective, call Emergency Services and give 5 Back Slaps.

3. If obstruction has not been cleared, follow by 5 Abdominal Thrusts.



unconscious, or if applied, u
both and nose p
000 send with hand
ACE • P
TIL BLOCKAGE
ot remove



Photoluminescent headers





