SHARSPGUARD®ECO

24L

Assembly Instructions



Match the lip of the lid and the lip of the base for correct assembly.



Put the temporary closure in place between use, press on mark for temporary closure.



To permanently close Press at back of door and slide thumbs forward following the arrows applying pressure to locking points.





Place hands either side of the lid and press until a "click" is heard.

Turn container 90 degrees and repeat, then press around the rim to ensure the lid is correctly attached to the base.



4 Lift or flick the central tab to open.



Do not forget to fill in the label following assembly and permanent closure.





REDEFINING SUSTAINABILITY



