


## Evolution LifeVac and Bleed Control Point



### Product Specifications

- **Life-Saving Device:** Combines a critical bleed kit with an advanced anti-choking device.
- **Regulated Medical Device** – CE marked, UKCA marked, MHRA regulated and MDR compliant
- **Comprehensive Emergency Coverage** – Together, these kits provide rapid response to two of the most time-critical incidents: catastrophic bleeding and airway obstruction, risks traditional first aid kits cannot adequately address.
- **Addresses Known High-Risk Scenarios** – Severe bleeding and choking are leading causes of preventable deaths. These solutions deliver targeted, life-saving intervention where standard kits may fall short.
- **Ideal** for any location seeking enhanced preparedness under increasing regulatory and safety expectations.

The **Evolution LifeVac and Bleed Control Point** is designed to address two potentially life-threatening emergencies: severe bleeding and airway obstruction. This combined provision supports prompt and effective first aid response, helping to enhance safety and emergency preparedness in high-risk environments. **Key Benefits**

- **Demonstrates compliance with evolving safety regulations** – Martyn Beadwell has highlighted a need for first aid supplies beyond those needed within general first aid. This first aid point addresses both critical bleeding and choking, two situations which could result in death if not dealt with promptly.
- **Proven life-saving equipment** – LifeVac has helped save more than 5,000 lives from choking incidents.
- **Enhances on-site emergency preparedness** – emergency situations can arise at any time in any location. So, having a variety of supplies available to treat them will enhance not only a workplace's preparedness but also build staff's confidence in crisis situations.
- **Includes a regulated medical device** – CE marked, UKCA marked, MHRA regulated and MDR compliant
- **Easily recognisable in emergency situations** – the boldly coloured Evolution cases are clearly labelled to ensure everyone is aware of their contents.
- **Supports current first aid provisions** – this first aid point can be added to pre-existing first aid measures to create an even more comprehensive cover to protect staff and visitors to the site.
- **Interlocking, quick release brackets** – wall-mountable, interlocking brackets not only allow your kits to be immediately visible and easily accessible but also allows customisation and reconfiguration to suit the space available.
- **Photoluminescent headers and bracket labels** – improve visibility and readiness by ensuring first aid points remain easy to locate in low-light conditions and by helping staff quickly identify and replace any missing kits.
- **Enables swift intervention during critical incidents** – the LifeVac allows rapid action in choking emergencies, while the critical injury kit provides essential supplies for life-threatening bleeds.
- **Promotes employee confidence and peace of mind** – knowing that both a critical bleed kit and life-saving choking devices are readily available helps staff feel secure, empowered, and prepared to respond effectively in serious, potentially life-threatening emergencies.
- **Cost-effective risk management solution** – this combination offers an affordable way to protect employees, visitors, and customers, while reducing potential legal claims and reputational damage.
- **Suitable for all UK workplaces** – ideal for schools, restaurants, transport hubs, warehouses, retail environments, and industrial sites where food is served, physical activity occurs, or injuries are possible, providing versatile safety coverage. **Why Have a LifeVac Anti-Choking Device on Site?** Choking can occur unexpectedly, affecting staff, visitors, or customers at any time. It is one of the UK's leading causes of accidental death, with brain damage possible within six minutes and death in under ten. A LifeVac device offers crucial protection, enabling a rapid and effective response before emergency services arrive. Keeping one on site, alongside a standard first aid kit, helps create a safer workplace for everyone. **Why Have a Bleed Control Kit on Site?** A

bleed control kit provides the essential tools, such as dressings, bandages, and tourniquets, needed to manage severe bleeding before emergency services arrive. Having a kit on site ensures staff can respond quickly and effectively, potentially saving lives and reducing the risk of long-term harm. Equipping a workplace with a bleed control kit demonstrates a commitment to employee safety, visitor protection, and overall workplace preparedness. **Kit**

**Contents** Bleed Control Kit:

- 1x Emergency Tourniquet
- 2x HypaStop Haemostatic Gauze Dressings, 10cm x 10cm
- 1x First Aid Guidance Leaflet
- 2x HypaTouch Nitrile Gloves (Pair)
- 2x HypaCover Trauma Dressings, Large
- 1x HypaGuard Foil Blanket
- 1x Clothing Cutters, 15cm

**LifeVac Kit**

- 1 x LifeVac User Guide Booklet
- 1 x LifeVac 'Instructions for Use' Leaflet
- 1 x LifeVac Unit.
- 1 x Paediatric Mask.
- 1 x Small Child/Adult Mask.
- 2 x Medium Adult (test) Mask.
- 1 x Large Adult Mask.
- 1 x A3 BLS Poster
- 1 x Window/Door Sticker

**Ideal For:**

- Catering and hospitality venues
- Schools, nurseries, and universities
- Factories, warehouses, and construction sites
- Public spaces and community centres
- Gyms, leisure centres, and sports facilities
- Transport hubs and vehicles



Official  
**LifeVac<sup>®</sup>**  
**Anti-Choking  
Device**

# First Aid Contents

## Bleed Control Kit Contents

Items	Qty
<b>Emergency Tourniquet</b>	<b>1</b>
<b>HypaStop Haemostatic Gauze Dressings, 10cm x 10cm</b>	<b>2</b>
<b>First Aid Guidance Leaflet</b>	<b>1</b>
<b>HypaTouch Nitrile Gloves (Pair)</b>	<b>2</b>
<b>HypaCover Trauma Dressings, Large</b>	<b>2</b>
<b>HypaGuard Foil Blanket</b>	<b>1</b>
<b>Clothing Cutters, 15cm</b>	<b>1</b>



## LifeVac® Choking Rescue Contents

Items	Qty
<b>Paediatric Mask</b>	<b>1</b>
<b>Small Child/Adult Mask</b>	<b>1</b>
<b>Medium Adult Mask</b>	<b>1</b>
<b>Large Adult Mask</b>	<b>1</b>
<b>LifeVac® Anti-Choking Device</b>	<b>1</b>
<b>Practice Mask (Medium Adult)</b>	<b>1</b>
<b>A3 BLS Poster</b>	<b>1</b>
<b>Window/Door Sticker</b>	<b>1</b>
<b>Instructions For Use Leaflet</b>	<b>1</b>
<b>LifeVac® Instruction Booklet</b>	<b>1</b>

**INSTRUCTIONS FOR USE IN THE CASE OF CHOKING**  
**If a casualty is a victim of choking**

- 1 Encourage the victim to cough.
- 2 If the cough is not effective, call Emergency Services and give 5 Back Slaps.
- 3 If obstruction has not been cleared, follow by 5 Abdominal Thrusts.



**Instructions For Use**

**1. Encourage the victim to cough.** If not effective call the emergency services and follow by 5 back slaps and 5 abdominal thrusts. Do not apply LIFEVAC until the obstruction has been cleared.

**2. If the cough is not effective, call Emergency Services and give 5 Back Slaps.**

**3. If obstruction has not been cleared, follow by 5 Abdominal Thrusts.**



unconscious, or if applied, u  
both and nose  
and with  
hand  
ACE • P  
TIL BLOCKAGE  
ot remove



# Photoluminescent headers





